IMPORTANT MEDICAL INFORMATION

Platinum/Iridium Upper Eyelid Implants for Lid Rehabilitation Post Facial Paralysis

Please read carefully prior clinical use!

GB

USA


Design	Product/Type	Reference Number	Material	Weight (g)	Intended Use	Description of the implant and the intended seating
	Upper Eyelid Implant Platinum/Iridium	4007 002- 4007 010	Platinum/ Iridium 90/10	0.6-2.2 (0.2g- increments)	KURZ Upper Eyelid Implants are intended for surgical implantation in the upper eyelid. They work by force of gravity to restore a functional blink mechanism of the eye.	The upper-eyelid platinum/iridium implant measures approx. 15-25,7 mm in length and 5.2 mm in width, depending on the version. The implants are 0,6-0,9 mm (REF 4007002) or 0.7-1.0 mm (REF 4007003-4007010) thin, and all edges are tapered and rounded in order to achieve good post-operative cosmetic results. The implant has a radius that is adjusted to the ocular bulb. Nine weights are available in total for customizing to individual requirements.

Important Medical Information KURZ- Platinum/Iridium Upper Eyelid Implants for Lid Rehabilitation due to Facial Paralysis (Bell's Palsy)

Indications:	Absence of or defective innervation of the orbicularis oculi muscle (lagophthalmos) due to peripheral facial paralysis, for example: • idiopathic facial paralysis lasting longer than 2 months; • concurrent with surgical reanimation of facial nerve, i.e. primary nerve suturing, rerouting, nerve interposition, VII/XII nerve anastomosis, free muscle graft; • injury to the facial nerve from trauma or tumor resection without possibility of nerve reanimation.
Contraindications:	• Known allergy to platinum/iridium • General wound-healing disorders • Lid closure unsuccessful with trial weight • Acute and chronic infectious diseases
Side Effects and Interactions:	Interaction with microwaves or magnetic resonance imaging (MRI); see Warnings.
Possible complications:	• Intolerability reactions • Migration • Corneal astigmatism • Tissue irritation • Contouring • Lid ptosis • Extrusion If a granuloma or other incompatibility reaction occurs, the necessary medical steps must be initiated immediately.
Surgical Technique:	(Non-binding recommendations, the surgeon and the anesthetist are responsible for the final decision regarding medication and dosage). Conjunctival and corneal anesthesia with 0.4% benoxinal solution. Placement of an eye covering cup. Local anesthesia of the skin and the subcutaneous tissue with 1% lidocaine hydrochloride and 1% adrenalin 100,000. Infiltration into the fold of the upper eyelid, followed by skin incision in the fold (approx. 3 mm longer than the implant). Severing of the subcutaneous tissue and the M. orbicularis oculi. It is crucial to work in a caudal direction when performing dissection. When working in a cranial direction, the M. levator palpebrae and its insertion in the tarsus are reached (injury hazard!). After sure identification of the tarsus, an appropriately large pocket is formed for accommodating the implant. The upper-eyelid implant is then fitted in place. The implant is fixed in place on the tarsus with 8-0 non-resorbable suture material (4 perforations). Suturing of the M. orbicularis oculi with 7-0 resorbable material. Closure of the skin incision with adhesive steri strips, removal in six days. No eye dressing required. Implantation can be permanent or temporary. In the latter case, the implant should be removed when the function of the buccal branch has been restored. In some cases, additional surgical measures, such as shortening of the lower eyelid, canthoplasty, eyebrow lifting or blepharoplasty, are required.
Implant Selection:	Preoperative determination of the respective implant (0.6 - 2.2 g) is possible by affixing a stainless steel weight specimen (REF 8000 111) by means of transparent adhesive plaster. Complete eyelid closure must be achieved. When the eyelids are opened, the upper eyelid should return to its original position. Caution: If too heavy an implant is selected, the palpebral fissure is diminished in size.
Sterility	KURZ Upper Eyelid implants are provided sterile. Packaging contents are sterile, provided that the packaging was not damaged or opened. KURZ implants are sterilized with gamma radiation in strictly controlled cycles.
Resterilization/ Reprocessing:	KURZ Upper Eyelid Implants are single-use products.
MRI:	MR-Conditional for 1,5 T; 3,0 T and 7,0 T for more detailed information about MRI please review www.kurzmed.com .
Warnings:	The physician must inform the patient about following: Patients with metallic implants shall not to be subjected to microwave radiation. In case of high levels of physical activeness, like sports, especially diving, the physician must inform the patient, that in the first few weeks after the surgery until up to 6 months, the eyelid has to be protected from external impacts, movement, pressure etc.
Storage:	Store at dry room temperature in unopened original package. Each implant bears a batch number and an expiration date and may not be used after that date.
Sale Restriction:	KURZ-prostheses are restricted to sale to or on the order of a specialist physician. Federal US law restricts KURZ Upper Eyelid Implants to sale to or on the order of a specialist physician.
Documentation:	The manufacturer recommends complete clinical, radiological and statistical documentation. Batch number, implant type and material must be entered into the patient's file, the surgical record or the implant card, using the adhesive labels provided with each implant.
 Manufacturer:	Heinz Kurz GmbH • Tuebinger Strasse 3 • 72144 Dusslingen • Germany • Tel.: +49 (0)7072 / 9179 – 0 • Fax: +49 (0)7072 / 9179 – 79 • Internet: www.kurzmed.com • E-Mail: info@kurzmed.de Address USA: Kurz Medical Inc. • 70 Chestnut Street • Shrewsbury, MA 01545 -USA • Phone: 508.841.5900 • Fax: 508.519.2672 • Toll Free: 866.449.8020 • Customer Service: orders@kurzmed.com
Version:	DMS 0001109_Rev.04-1 –9501430

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EN

SYMBOLS

	Manufacturer		Item number
	Expiration date		Batch code
	Sterilized using irradiation		Do not use if package is damaged
	Keep away from sunlight		Do not resterilize
	Do not reuse!		Single sterile barrier system
	Store in a dry place		Single sterile barrier system with protective packaging inside
	Consult instructions for use		Single sterile barrier system with protective packaging outside
	Fragile, handle with care		Medical device
			Caution