GENERAL INFORMATION FOR STAPEDIOPLASTY PROSTHESES AND ACCESSORIES

en	GENERAL INFORMATION FOR STAPEDIOPLASTY PROSTHESES AND ACCESSORIES	4
de	ALLGEMEINE HINWEISE ZU PROTHESEN FÜR DIE STAPESPLASTIK UND ZUBEHÖR	7
fr	RECOMMANDATIONS GÉNÉRALES CONCERNANT LES PROTHÈSES DE STAPÉDOPLASTIE ET LEURS ACCESSOIRES	10
it	AVVERTENZE GENERALI SULLE PROTESI PER LA STAPEDIOPLASTICA E RELATIVI ACCESSORI	13
es	INDICACIONES GENERALES RELATIVAS A LAS PRÓTESIS PARA ESTAPEDIOPLASTIA Y ACCESORIOS	16
pt	INDICAÇÕES GERAIS RELATIVAS A PRÓTESES PARA ESTAPEDOPLASTIA E ACESSÓRIOS 1	19
ru	ОБЩИЕ УКАЗАНИЯ ПО ПРОТЕЗАМ ДЛЯ СТАПЕДОПЛАСТИКИ И ПРИНАДЛЕЖНОСТЯМ	22
cs	VŠEOBECNÉ POKYNY K PROTÉZÁM PRO PLASTIKU TŘMÍNKU A PŘÍSLUŠENSTVÍ	25
sr	OPŠTE INFORMACIJE O PROTEZAMA I PRIBORU ZA STAPEDOPLASTIKU	28
pl	OGÓLNE WSKAZÓWKI DOTYCZĄCE PROTEZ DO STAPEDOPLASTYKI I WYPOSAŻENIA DODATKOWEGO	31
da	GENERELLE HENVISNINGER TIL PROTESERNE TIL STAPESPLASTIK OG TILBEHØR	34
sv	ALLMÄNNA ANVISNINGAR GÄLLANDE PROTESER FÖR STAPESPLASTIK OCH TILLBEHÖR	37

no	GENERELLE MERKNADER OM PROTESER FOR PLASTISK KIRURGI AV STIGBØYLE OG TILBEHØR	40
lv	VISPĀRĒJAS NORĀDES PAR KĀPSLĪŠA PLASTIKAS PROTĒZĒM UN PIEDERUMIEM	43
lt	BENDROJI INFORMACIJA APIE STAPEDOPLASTIKAI SKIRTUS PROTEZUS IR PRIEDUS	46
nl	ALGEMENE OPMERKINGEN OVER PROTHESES VOOR DE STAPESPLASTIEK EN ACCESSOIRES	49
tr	STAPEDİOPLASTİ PROTEZLERİ VE AKSESUARLARI İÇİN GENEL BİLGİLER	52
el	ΓΕΝΙΚΕΣ ΥΠΟΔΕΙΞΕΙΣ ΣΧΕΤΙΚΑ ΜΕ ΤΙΣ ΠΡΟΘΕΣΕΙΣ ΓΙΑ ΤΗΝ ΠΛΑΣΤΙΚΗ ΑΝΑΒΟΛΕΑ ΚΑΙ ΤΑ ΠΑΡΕΛΚΟΜΕΝΑ	55
sl	SPLOŠNA NAVODILA ZA PROTEZE ZA STAPEDOTOMIJO IN DODATNA OPREMA	58
hr	OPĆE INFORMACIJE O PROTEZAMA I PRIBORU ZA STAPEDOPLASTIKU	61
sk	VŠEOBECNÉ POKYNY K PROTÉZAM PRE STAPEDOPLASTIKU A PRÍSLUŠENSTVO	64
ro	INDICAȚII GENERALE PRIVIND PROTEZELE PENTRU STAPEDECTOMIE ȘI ACCESORIILE ACESTORA	67
zh	关于镗骨造型术用假体和配件的一般提示	70
ja	アブミ骨形成術用プロテーゼおよび付属品に関する一般情報	73

Please read carefully before clinical use!

INTENDED USE: In principle, KURZ middle ear prostheses provide for partial or total replacement of the human auditory ossicular chain when surgically indicated. The objective is to transmit sound energy along the mechanical path from the eardrum to the oval window of the cochlea with the least possible losses.

Accessories

KURZ recommends the use of original accessories such as

- SteadyCrimp Forceps REF 8000188
- Soft-Clip Hook REF 8000127
- KURZ Meter with Tray REF 8000106

Indications

The following underlying diseases may lead to an indication of a passive middle-ear prosthesis for reconstruction of the human auditory ossicular chain:

- Otosclerosis / congenital stapes fixation
- Malformation of the middle ear with fixed foot plate
- Revision surgery due to inadequate hearing improvement, e. g. due to dislocation of a (previously implanted) prosthesis

Contraindications

- Allergy to titanium
- Situations where conservative treatment approaches are adequate
- Acute middle ear inflammation that can result in dislocation of the prosthesis
- Acute and chronic infectious diseases
- General wound healing impairments
- Important! Any opening-up of the inner ear (such as perforation of the footplate of the stapes, partial or complete removal of the stapes footplate (stapedioplasty) is contraindicated in cases of inflammation in the acoustic meatus and/or the middle ear due to the risk of the inflammatory process spreading into the inner ear
- Concerns only the product NITIFLEX stapes prosthesis made of nickel-titanium alloy (Nitinol): Allergy to nickel or nickel-titanium alloys (nitinol). A preoperative allergy test is recommended to minimize the risk of an allergic reaction! In the case of patients with known hypersensitivity to nickel or nickel-titanium alloys (nitinol), the use of an implant containing any of these materials must be avoided.

Potential complications / Undesirable side effects

Possible adverse events / injury may occur during or after surgery. Fine bony structures are touched and mobilized as part of the surgical procedure and this may result in surgically induced trauma or infection. Such damage may be irreversible or correctable only through revision surgery.

- Postoperative dislocation of the prosthesis
- Incus necrosis
- Recurrent or postoperative middle-ear inflammation
- Vertigo
- Labyrinthine hearing loss following intraoperative trauma
- Tissue irritation, scar formation, granuloma
- Perilymphatic fistula
- Perforation of the tympanic membrane
- Inner ear damage up to and including deafness (Surditas)
- Tinnitus
- Irritation or even damage to the facial nerve up to and including paresis of the facial nerve (e. g. intraoperative injury of the corda tympani)
- Incus subluxation
- Floating footplate
- Lack of hearing improvement

Intended User

KURZ middle ear prostheses and accessories are restricted to be sold to and used by specialized ENT Surgeons.

PRECAUTIONS / WARNINGS

- Intraoperative

Unintentional bending of the implant must be absolutely avoided when removing it from the primary packaging, since this may cause defective functioning of the implant itself. Extreme care must be taken in the selection of the prosthesis length for the implantation to prevent later problems such as implant dislocation or the occurrence of inner ear symptoms, such as vertigo, for example.

- Postoperative

Should a granuloma or perilymphatic fistula develop after surgery, appropriate medical procedures must be initiated immediately.

- Surgical Technique

KURZ middle ear prostheses and accessories are meant to be used in an operating room.

The surgical procedure (exemplary for Stapedotomy) can be described in short form as follows:

- Stapedioplasty (under microscopic view) is preferably conducted under general anesthesia
- The middle ear is opened endaurally or retroauricularly
- Formation of a tympano-meatal flap
- Releasing of the annulus tympanicus
- Folding back of the tympanic membrane, avoiding damaging of the chorda tympani and proceeding into the tympanic cavity
- For clear view of the pyramid the postero-superior canal wall is curetted out
- Palpative examination of the ossicles and footplate including palpation of the stapes capitulum to confirm the diagnosis
- Separation of the incudo-stapedial joint plus tenotomy
- Partial or complete crurotomy and removal of the stapes capitulum
- Fenestration of the footplate (Stapedotomy) preferably by laser
- The specialized ENT surgeon chooses the prosthesis in accordance with the anatomical condition and ossicular defect
- Insertion of the stapes prosthesis into the fenestrum and coupling the head of the prosthesis to its final position with the long process of incus by crimping, clipping
- Confirming the appropriate fixation of the prosthesis to the incus
- Sealing the piston with tissue around the fenestration
- Repositioning of the tympano-meatal flap and tamponating the ear canal

- MRI Safety Information

MR Conditional for 1.5 Tesla, 3.0 Tesla and 7.0 Tesla. For more detailed information, see Internet: www.kurzmed.com

- Warning

The patient must be informed by the physician about the following aspects. It is not permissible to expose patients with metal implants to microwave radiation. Strong fluctuations of the ambient pressure (immersion, diving headfirst into water, explosions, fireworks, etc.), sport (e. g. types of contact sports) as well as sauna visits (during the healing phase) may damage middle-ear structures and/or lead to disturbances of the sense of hearing and equilibrium, and are therefore to be avoided.

- Sterility

KURZ middle ear prostheses are delivered in sterile condition. The contents of the packaging are sterile, provided the packaging has not been damaged or opened. KURZ prostheses are sterilized with gamma radi-

tion in strictly controlled cycles. The storage packaging is permitted to be opened only immediately before the operation.

- Resterilization / Reprocessing

The prosthesis is intended for one-time use only. Reprocessing / resterilization is not permitted.

For resterilization of accessories the instructions for "Reprocessing Instructions for Cleaning, Disinfection and Sterilization of KURZ Instruments" must be followed.

The end of accessories product lifetime is significantly determined by the wear, the treatment process, the chemicals used and any possible damage caused by use. If disregarded, any liability is excluded.

- Storage

Storage conditions for unopened original packaging: Store in a dry place at room temperature, protected from direct sunlight. Brief fluctuations of temperature and humidity are permissible. Every prosthesis is marked with a LOT number and an expiration date. Once the expiration date has passed, it is no longer permitted to be implanted.

- Disposal

Please observe the country-specific regulations on waste disposal.

- Documentation

The manufacturer recommends that the LOT-number and the prosthesis type be recorded on the medical record, the surgical report and the implant pass by means of the adhesive labels.

Sales Restriction

US Federal Law restricts KURZ prostheses to sale by or on the order of a physician.

Incidents

All serious incidents that have occurred in connection with the product must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Summary of Safety and Clinical Performance

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en

SYMBOLS

	Manufacturer		Item number
	Expiration date		Batch code
	Sterilized using irradiation		Do not use if package is damaged
	Keep away from sunlight		Do not resterilize
	Do not reuse!		Single sterile barrier system
	Store in a dry place		Single sterile barrier system with protective packaging inside
	Consult instructions for use		Single sterile barrier system with protective packaging outside
	Fragile, handle with care		Medical device
			Caution

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