SteadyCrimP Forceps (REF 8000188)

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0004698_Rev05_9501635_ 2024-08_de_en_fr_es

en **INSTRUCTIONS FOR USE FOR STEADYCRIMP FORCEPS**
Read carefully prior to clinical use!

INTENDED PURPOSE: The SteadyCrimp Forceps are a passive, reusable device which is used intraoperatively and surgically invasive to crimp / fix a KURZ crimpable Stapes prosthesis onto the long process of the incus or the manubrium of the malleus.

Application / Description

The SteadyCrimp Forceps are designed to attach a crimp-stapes prosthesis onto the long process of the incus or the manubrium of the malleus during a stapesplasty.

The SteadyCrimp Forceps are a reusable medical device and can be reprocessed.

In the following we would like to provide you with some important information:

- Handle the instrument with care
- Avoid impact, bending and other mechanical overloading
- Do not use lubricants or surface oils

Intended User: Specialized ENT surgeons

Intended Patient Population

The intended patient population, patient selection criteria correspond to the intended patient population, patient selection criteria of a KURZ Stapedioplasty prosthesis. All patients who require a Kurz middle ear implant are suitable for the procedure after a careful differential diagnosis. Children must have reached a certain level of maturity. The degree of maturity is decided by the physician performing the procedure. According to the clinical data, there is a very successful result in children from the age of 5 years.

Indications

The indication for use corresponds to the indication for use of a crimpable KURZ Stapedioplasty prosthesis.

Contraindications

The contraindications correspond to the contraindications of a crimpable KURZ Stapedioplasty prosthesis. No further contraindications.

Potential complications / Undesirable side effects: None.

PRECAUTIONS / WARNINGS

- If crimping is too tight, this can lead to necrosis of the incus
- If crimping is too loose, this can lead to dislocation of the prosthesis or to an insufficient audible outcome
- Be careful when removing the instrument to avoid damage
- Clean, disinfect and sterilize the instrument completely before initial use

Disposal

Please observe the country-specific regulations on waste disposal.

Reprocessing Instructions for the Cleaning, Disinfection and Sterilization of KURZ Instruments

For detailed information on reprocessing, please refer to the "Reprocessing instructions for instruments" in the respective national language.

Incidents

All serious incidents that have occurred in connection with the product must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Summary of Safety and Clinical Performance


https://ec.europa.eu/info/index_en

Technical Manager

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
DISASSEMBLY / ASSEMBLY / CLEANING


DISASSEMBLY:


a) Assembled instrument.


b) Open the lock by turning the latch. **Info:** The latch slides on a beveled surface of the pin and is thus lifted, releasing the lock. **Attention:** Do not lift the latch with force.

c) Carefully remove the thumb grip, including the upper part, from the lower part in the direction of the arrow.


d) Detach the thumb grip from the upper part by moving it sideways.


e) **Info:** The jaw section with joint cannot be disassembled further.


ASSEMBLY:


a) Join the thumb grip to the upper part by moving it sideways.

b) Turn the oblong slot of the latch to align it with the slot of the lower part.


c) Press the thumb grip with the upper part in the direction of the arrow and allow to engage into the lower part.

d) Check the locked thumb grip with top part for:
• Secure locking, by pulling the thumb grip upwards in the direction of the arrow.
• Functioning of the pliers by opening and closing the jaws of the pliers.

CLEANING:

Disassembled instrument and latch position during cleaning.

SYMBOLS

	Manufacturer		Consult instructions for use		Single sterile barrier system
	Expiration date		Fragile, handle with care		Single sterile barrier system with protective packaging inside
	Sterilized using irradiation		Item number		Single sterile barrier system with protective packaging outside
	Keep away from sunlight		Batch code		Medical device
	Do not reuse!		Do not use if package is damaged		Caution
	Store in a dry place		Do not resterilize		

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