

Tympanoplasty Prostheses



MRP Malleus Replacement



HEINZ KURZ GMBH
TUEBINGER STR. 3
72144 DUSSLINGEN
GERMANY

Table of Contents

1 About this Document	3		
1.1 Symbols Glossary	3		
1.2 Safety Information Marking.....	3		
1.3 Additional Information	4		
1.4 Safety-related Changes.....	4		
2 Important Safety Information	4		
3 Product Codes / REF	4		
4 Scope of Delivery	4		
5 Packaging and Sterility	4		
6 Product Description	5		
6.1 General Information.....	5		
6.2 Structure and Operation.....	5		
6.3 Materials with Potential Patient Contact.....	5		
6.4 Accessories	5		
6.5 Other Devices to be Used in Combination with the Device	5		
7 Intended Use	5		
7.1 Intended Purpose.....	5		
7.2 Indications.....	5		
7.3 Contraindications	5		
		7.4 Patient Target Group	6
		7.5 Intended User	6
		7.6 Expected Lifetime	6
		7.7 Intended Place of Use	6
		8 Expected Clinical Benefit	6
		9 Possible Complications and Side Effects	6
		10 Combining with Other Procedures	6
		11 Shelf Life and Storage.....	6
		12 Processing.....	7
		13 Application Instructions	7
		13.1 Required Equipment and Materials.....	7
		13.2 Preparation of the Patient	7
		13.3 Preparing the Prosthesis.....	7
		13.4 Placing the Prosthesis.....	8
		13.5 Removing the Prosthesis	8
		14 Aftercare	8
		15 Instructing the Patient	9
		16 Implant Card	9
		17 Disposal	9
		18 Warranty	9
		19 Specifications	10

1 About this Document

1.1 Symbols Glossary

Symbol	Description
	Caution: Consult Instructions for Use
	Caution!
	Fragile; handle with care
	Do not use if package is damaged
	Keep away from direct sunlight
	Keep dry
	Use-by date
	Sterilized using irradiation
	Do not reuse
	Do not resterilize
	Single sterile barrier system with protective packaging inside
	MR conditional
	Medical device
	Catalog number
	Batch code
	Unique Device Identification (UDI)
	Quantity per packaging unit
	Manufacturer
	Date of manufacture
	(USA) Caution: Federal Law restricts this device to sale by or on the order of a physician.
	Consult Instructions for Use. The Instructions for Use are provided in electronic form (e-labelling).
	Patient name
	Date of implantation
	Name of the implanting healthcare institution / provider
	Patient information website

Table 1: Symbols Glossary

1.2 Safety Information Marking

WARNING

Non-compliance may result in serious injuries, serious deterioration of the general condition or the death of the patient, user, or a third party.

NOTICE

Product damage or other damage may occur in case of non-compliance.

1.3 Additional Information

This document is available in electronic form on the manufacturer's website. If required, a printed copy of this document can be requested from the manufacturer.

Download link for these Instructions for Use: ¹⁾	www.kurzmed.com/en/ifu/tym7.html
 Download link for the Patient Information Document: ¹⁾	www.kurzmed.com/en/pi/tym.html
Disclaimer for the availability of the SSCP	As a general rule: The SSCP will only be made available after the product has been authorised in accordance with REGULATION (EU) 2017/745 (MDR). The implementation described here does not apply until the corresponding module of the Eudamed database comes into force. Until then, the SSCP is available at the following download link: www.kurzmed.com/en/sscp/tym.html
Summary of Safety and Clinical Performance (SSCP): ¹⁾	https://ec.europa.eu/tools/eudamed To search for the product-specific SSCP, enter the basic UDI-DI of the product.
Basic UDI-DI (device identifier):	++EHKM0017D
International addresses:	https://www.kurzmed.com/en/contact.html

¹⁾ Updated on an ongoing basis. Other language versions are also available there.

For the complete UDI (UDI-PI) please refer to the product label.

1.4 Safety-related Changes

Document number	Edition date	Safety-related Changes
0005959_01	2024-10	Complete revision
0005959_02	2026-02	Product discontinued: Malleus Handle Cavity Bending Pliers

2 Important Safety Information

WARNING

- Before using the product: Read the Instructions for Use for the product and for all products used in combination. Follow and save the Instructions for Use.
Otherwise there are risks for the health of your patient.
- Do not disassemble or modify the product.
Otherwise there are risks to the health of your patient.

ATTENTION: In case that any serious incident has occurred in relation to the device the incident should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

3 Product Codes / REF

[▶ Specifications, page 10]

4 Scope of Delivery


MRP Malleus Replacement (Tympanoplasty Prosthesis)	1 x prosthesis 1 x implant card 4 x product label
---	---

5 Packaging and Sterility

MRP Malleus Replacement (Tympanoplasty Prosthesis)	The product is sterile (sterilized by radiation). Packaging: Single sterile barrier system with protective packaging inside (prosthesis in plastic triangular box and hard blister) + outer packaging (folding box)
---	--

6 Product Description

6.1 General Information



- A Replacement malleus handle for linking with a KURZ partial/total prosthesis
- B Y-shaped shaft
- C Pins for anchoring in the auditory canal wall

Illustration 1: MRP Malleus Replacement

[▶ Specifications, page 10]

6.2 Structure and Operation

MRP Malleus Replacement (Tympanoplasty Prosthesis)	Prostheses which are inserted to partially or completely replace middle ear structures involved in sound conduction.
---	--

6.3 Materials with Potential Patient Contact

The following table lists all implant materials that the user or patient may come into contact with during application.

Product (part)	Material	Contact person
MRP Malleus Replacement (Tympanoplasty Prosthesis)	100% titanium	Patient

Not made with natural rubber (latex).

No products made with natural rubber (latex) are used in the production process.

ATTENTION: Do not use the product if the patient has known intolerances / allergies to the materials used.

6.4 Accessories

Accessories (separate instructions for use):

- KURZ Precise Cartilage Knife Set (REF 8000 155)
- Cartilage Forceps Schimanski Design (REF 8000 193)

6.5 Other Devices to be Used in Combination with the Device

MRP Malleus Replacement is intended for use together with various KURZ partial/total prostheses.

Compatibility: [▶ Specifications, page 10]

7 Intended Use

7.1 Intended Purpose

MRP Malleus Replacement (Tympanoplasty Prosthesis)	KURZ middle ear prostheses are intended for the partial or total surgical replacement of the ossicular chain of the human middle ear. The objective is the restoration of the mechanical transfer of sound from the tympanic membrane to the oval window of the cochlear with the least impairment of hearing.
---	---

7.2 Indications

- Chronic otitis media with functional impairment of the ossicular chain
- Traumatic injury to the ossicular chain
- Congenital malformations of the middle ear
- Revision surgery due to inadequate hearing improvement (e.g., due to dislocation of a previously implanted prosthesis)

7.3 Contraindications

- Known sensitivity or allergy to titanium
- Complications or sequelae of unresolved otitis media, such as intracranial abscess, meningitis, lateral sinus thrombosis, malignancies, or patient-specific systemic disease
- Acute middle ear inflammation

- Impaired wound healing

7.4 Patient Target Group

The product is suitable for use in the following patient groups:

- Children and youth
- Adults
- Patients of any gender

7.5 Intended User

The intended user is a physician with experience in treating similar cases with this product or with comparable products or a physician with the following specialty:

- ENT (otorhinolaryngology)

7.6 Expected Lifetime

No product-specific restrictions.

7.7 Intended Place of Use

- Operating theatre

It is the responsibility of the user to decide on a case-by-case basis which precautions must be taken for any complications that may arise.

8 Expected Clinical Benefit

According to the clinical evaluation, the product can be used safely and effectively for treatment according to the indications mentioned.

9 Possible Complications and Side Effects

- Implant migration
- Implant extrusion
- Lateralisation of the implant
- Sensorineural hearing loss
- Infection
- Dizziness
- Periprosthetic fibroses
- Periprosthetic cholesteatom formation

10 Combining with Other Procedures

WARNING

- Laser therapy, argon plasma coagulation, high-frequency surgery, and other procedures, the effect of which is due to heat: Do not use those methods directly on the product. Otherwise, injury to the tissue and product damage are possible.
- Apart from the specific provisions on MRI safety, the following applies: Do not expose the product to diagnostic or therapeutic electromagnetic radiation. Otherwise there are risks to the health of your patient.
- The product is MRI conditional. Use the product in MRI fields only as per specification. Possible consequences of using the product in MRI fields outside the specifications include: Heating of the product, electrostatic discharges, consequential damages caused by the application of force to the product, errors in the imaging (also in the surrounding tissue)

For important information about MRI see:

<http://www.kurzmed.com/de/mr-information.html>

11 Shelf Life and Storage

For date of expiry, see the product label.

Store the product in unopened original packaging.

Store the product in a dry place and protect it from sunlight.

12 Processing

⚠ WARNING

- Single use product: Do not process (e.g., clean, disinfect, sterilize), resterilize or reuse the product.
This is the only way to ensure the product is germ-free and functional. Due to the mechanical properties of the product, processing or resterilization could lead to material degradation.

13 Application Instructions

⚠ WARNING

- Do not use the product if the packaging or the product is damaged or expired.
This is the only way to ensure the product is germ-free and functional.
- Only remove the product from storage packaging immediately before use. When the product is removed from the packaging, observe the relevant hygienic regulations.
Otherwise there are risks to the health of your patient.

NOTICE

- Always grasp, transport and manipulate the prosthesis with a suitable suction device or with appropriate forceps or tweezers. Ensure that the prosthesis shaft is not inadvertently deformed or the prosthesis is not damaged in any other way.
Otherwise the function of the prosthesis may be impaired.

Ensure the presence of hygienic / sterile conditions needed for the intervention.

It is placed as part of a type III tympanoplasty (ossicular reconstruction).

Perform the intervention under appropriate visual supervision.

ATTENTION: Also observe the instructions for use of the KURZ partial/total prosthesis used.

13.1 Required Equipment and Materials

As usual for a type III tympanoplasty.

- Compatible KURZ partial/total prosthesis [▶ Specifications, page 10]

The manufacturer recommends using the following products:

- KURZ Precise Cartilage Knife Set (REF 8000 155)
- Cartilage Forceps Schimanski Design (REF 8000 193)

13.2 Preparation of the Patient

As usual for a type III tympanoplasty.


Endaural or retroauricular access to the middle ear.

13.3 Preparing the Prosthesis




1. Open the sterile packaging.
2. Apply drops of sterile saline solution on the openings of the protective packaging. In this process, ensure that the perforations in the lid are also coated in saline solution so that liquid can penetrate the protective packaging.
3. Carefully remove the prosthesis from the protective packaging. ATTENTION: Do not grasp the prosthesis by the shaft to avoid bending the prosthesis.


13.4 Placing the Prosthesis




1. Drill two holes into the auditory canal wall. Choose the position of the holes according to the overall anatomical situation/condition of the auditory canal wall.
Hole diameter: 0.6 mm
Depth: Approx. 2 mm
Centre-to-centre distance: 3.0 mm
ATTENTION: During drilling, flush with water for cooling purposes.




2. Insert the two pins of the Malleus Replacement Prosthesis into the holes.



3. Adapt the Malleus Replacement Prosthesis to the anatomical conditions. For this purpose, carefully bend the shaft of the Malleus Replacement Prosthesis. Then insert the KURZ partial/total prosthesis. See the instructions for use of the partial/total prosthesis.



4. Stabilise the partial/total prosthesis with the Malleus Replacement Prosthesis. For this purpose, position the replacement malleus handle of the Malleus Replacement Prosthesis in the malleus handle cavity of the partial/total prosthesis' head plate.



5. Completely cover the Malleus Replacement Prosthesis and the partial/total prosthesis' head plate against the tympanic membrane with a graft (cartilage disc, thickness approx. 0.3 - 0.5 mm).

13.5 Removing the Prosthesis

The prosthesis is intended to remain in the body. However, should it nevertheless be necessary to remove the prosthesis: Before removing the prosthesis: Loosen any adhesions.
Follow-up treatment at the discretion of the treating physician.

14 Aftercare

- Follow-ups as indicated by the treating physician.

15 Instructing the Patient

The instruction to the patient must include:

⚠ WARNING

- Protect the auditory canal from water penetration.
Otherwise there is a risk of inflammation / infection of the middle ear.
- Avoid severe fluctuations in ambient pressure (e.g. diving, jumping head first into water, explosions).
Failure to do so may result in injury to the tympanic membrane/ossicles, which can lead to hearing and balance disorders.


ATTENTION: Also inform the patient about the consequences of combining with other procedures.

[▶ Combining with Other Procedures, page 6]

Implant card: [▶ Implant Card, page 9]

16 Implant Card

ATTENTION: Complete the implant card before the patient is discharged from hospital and give it to the patient.



1. Stick one of the product labels provided into the designated box on the implant card. Complete all other sections.

The implant card must be presented at every radiological examination.

17 Disposal

⚠ WARNING

- The product was in contact with potentially infectious substances of human origin. Clean/pack the product for disposal according to the specific contamination risk. Dispose of the product in accordance with procedures for hazardous waste in hospitals.
Otherwise there is a risk of infection for the user and for third parties.

Disposal must be in accordance with national disposal regulations and pursuant to the corresponding risk class.

18 Warranty

The reliability of the product's material and design at the time of shipment is guaranteed. The manufacturer does not know either the diagnosis of the patient or the nature of the application and has no influence on the conditions under which the product is used. The storage conditions after delivery of the product are also beyond the manufacturer's area of responsibility.

Due to biological and individual differences, no product is 100% effective under all circumstances.

Therefore, the manufacturer cannot guarantee a positive effect or the absence of negative effects for product application. The medical staff must use the product on the basis of their medical training and experience, and they are responsible for correct application.

The warranty (repair or replacement) applies only if the product is used in accordance with these Instructions for Use (for instruments, particularly with regard to handling, cleaning, sterilization and maintenance); the warranty period starts on the delivery date.

If you have reason to believe that a new product is faulty, please contact the Customer Service in writing immediately and provide as detailed a description as possible of the fault, the REF (product code), and the LOT (batch code) and/or series number. All allegedly defective products must be returned to us for inspection. Instruments have to be completely cleaned and sterilized, appropriate documentation must be enclosed with the return.

If the manufacturer finds that despite all due care the product was defective at the time of delivery, he will repair the product or replace it promptly. If repair or replacement of the product is not possible, the buyer has the right to cancel the purchase or to reduce the payment, but by a maximum of the purchase price amount.


Additional claims or those not mentioned here due to defect, and other claims regardless of the legal reason, including those based on illegal acts and for compensation of immaterial damages against the manufacturer, his agents, dealers and suppliers, are excluded unless existing law is contrary to the liability exclusion, e.g. in cases of intent or gross negligence or in the event of physical injury.

All claims based on the consequences of non-compliance with the Instructions for Use, including specified indications, contraindications, warnings, instructions, application, storage and off-label use, as well as the consequences of a combination with third-party products are excluded.

Furthermore, all claims that result from the use of products that have expired, or were used despite the obvious damage to the packaging, or resterilized and/or recycled contrary to the Instructions for Use, are excluded.

No one is allowed to change the above conditions, make further warranty or liability declarations, or guarantee any properties that surpass those specified in the Instructions.

19 Specifications

	Name	REF	Compatible KURZ partial/total prostheses
	MRP Malleus Replacement	1006 960	• MNP Malleus Notch Partial • MNP Malleus Notch Total