(DE)	ALLGEMEINE HINWEISE ZU PROTHESEN FÜR DIE OSSIKELPLASTIK.....	2
(EN)	GENERAL INFORMATION FOR OSSICULOPLASTY PROTHESES.....	4
(FR)	INFORMATIONS GENERALES SUR LES PROTHESES POUR OSSICULOPLASTIE.....	6
(IT)	INFORMAZIONI GENERALI SULLE PROTESI PER LA OSSICULOPLASTICA.....	8
(ES)	OBSERVACIONES GENERALES SOBRE LAS PRÓTESIS PARA OSSICULOPLASTIA.....	10
(PT)	INDICAÇÕES GERAIS RELATIVAS A PRÓTESES PARA OSSICULOPLASTIA.....	12
(RU)	Общие указания по работе с протезами для оссикулопластики.....	14
(CZ)	VŠEOBECNÉ POKYNY K PROTÉZÁM PRO OSIKULOPLASTIKU.....	16
(SRB)	OPŠTE NAPOMENE O PROTEZAMA ZA REKONSTRUKCIJU SLUŠNIH KOŠČICA.....	18
(PL)	OGÓLNE WSKAZÓWKI DOTYCZĄCE PROTEZ DO PLASTYKI KOSTECZEK SŁUCHOWYCH.....	20
(IS)	ALMENNAR RÁÐLEGGINGAR UM GERVILIDI FYRIR HEYRNABEINSPLASTID.....	22
(DK)	GENERELLE ANVISNINGER VEDRØRENDE PROTESER TIL OSSIKULOPLASTIK.....	24
(FI)	OSSIKKELIPLASTIIKASSA KÄYTETTYJÄ PROTEESEJÄ KOSKEVIA YLEISIÄ OHJEITA.....	26
(SE)	ALLMÄNNA ANVISNINGAR FÖR PROTESER FÖR OSSIKULOPLASTIK.....	28
(NO)	GENERELLE MERKNADER OM PROTESER FOR OSSIKELPLASTIKK.....	30
(TR)	OSSİKELPLASTİK İÇİN PROTEZLERE İLİŞKİN GENEL BİLGİLER.....	32
(NL)	ALGEMENE AANWIJZINGEN VOOR PROTHESEN VOOR DE OSSICULUMPLASTIEK.....	34
(GR)	ΓΕΝΙΚΕΣ ΥΠΟΔΕΙΞΕΙΣ ΓΙΑ ΤΗΝ ΠΛΑΣΤΙΚΗ ΑΚΟΥΣΤΙΚΩΝ ΟΣΤΑΡΙΩΝ.....	36
(HU)	ÁLTALÁNOS TANÁCSOK A KIS CSONTOK PROTETIZÁLÁSÁHOZ.....	38
(SI)	SPLOŠNI NAPOTKI ZA PROTEZE ZA KOŠČIC SREDNJEGA UŠESA.....	40
(RO)	INDICAȚII GENERALE PRIVIND PROTEZELE PENTRU PROTEZAREA OSICULELOR.....	42
(EE)	ÜLDISED NÕUANDED KUULMELUUKESTE PLASTIKAKS KASUTATAVATE PROTEESIDE KOHTA.....	44
(LT)	BENDROSISIOS PASTABOS DĖL OSIKULOPLASTIKOS PROTEZŲ.....	46
(CN)	听小骨假体概况.....	48

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Read carefully prior to clinical use!

Intended Use

KURZ middle ear prostheses are intended for the partial or total surgical replacement of the ossicular chain of the human middle ear. The objective is the restoration of the mechanical transfer of sound from the tympanic membrane to the oval window of the cochlear with the least impairment of hearing.

Indications

Implantation of a passive middle ear prosthesis to reconstruct the human ossicular chain can be indicated on account of the following underlying diseases:

- Chronic middle ear inflammation (also following removal of a tumour, e.g. cholesteatoma) with functional impairment of the ossicular chain
- Otosclerosis / congenital stapes fixation
- Traumatic injury to the ossicular chain
- Malformation of the middle ear
- Revision surgery due to inadequate hearing improvement, e.g. due to dislocation of a (previously implanted) prosthesis

Contraindications

- Allergy to titanium
- Cases that may be adequately treated by well-established, conservative methods
- Acute middle ear inflammation that could cause prosthesis dislocation
- Acute and chronic infectious diseases
- General impairment of wound healing

Important! Perforation of the footplate of the stapes and thus opening-up of the inner ear (stapedioplasty) is contraindicated in cases of inflammation in the acoustic meatus and/or the middle ear due to the risk of the inflammatory process spreading into the inner ear.

Only for nickel-titanium-prostheses (nitinol)!

- Allergy to nickel or nickel-titanium alloys (nitinol)
A preoperative allergy test is recommended to minimize the risk of an allergic reaction. In patients with a known allergy to nickel or nickel-titanium alloys (nitinol), do NOT use an implant containing these materials.

Possible Complications / Unanticipated adverse effects

Adverse effects or injuries may occur during or after the operation. During a surgical procedure delicate bone structures are touched and moved which may result in operative trauma or infection. These injuries may be irreversible or only reversible with revision surgery:

- Postoperative dislocation of the prosthesis
- Incus necrosis
- Recurrent middle ear inflammation
- Vertigo
- Labyrinthine hearing loss due to intraoperative trauma
- Tissue irritation, development of scar tissue, granuloma
- Perilymphatic fistula
- Perforation of the tympanic membrane
- Various degrees of inner ear injury up to deafness
- Tinnitus
- Irritation or even damage to the facial nerve, including facial paralysis

Precautions / Warnings

Intraoperative

Accidental bending of the implant during removal from its primary package must be avoided at all costs in order to prevent functional damage. Extreme care must

be exercised when choosing the proper length of the prosthesis to be implanted to avert postoperative problems such as implant dislocation or inner ear symptoms such as vertigo. In tympanoplasty, implants that touch the tympanic membrane must be covered with a cartilage slice (0.3-0.5 mm).

Postoperative

Should a granuloma or a perilymphatic fistula develop after surgery, appropriate medical procedures must be initiated immediately.

MRI Safety Information

MR Conditional at 1.5 T; 3.0 T and 7.0 T, for detailed Information please refer to the MR Information on www.kurzmed.com.

Warning

The physician must inform the patient about the following: Patients with metallic implants are not to be subjected to microwave radiation. Severe variations in ambient pressure (scuba diving, diving head-first, explosions, fireworks, etc.) can damage the middle ear structures and/or cause hearing loss and/or equilibrium disorder and should be avoided.

Sterility

KURZ middle ear prostheses are provided sterile. Packaging contents are sterile, provided that the packaging was not damaged or opened. KURZ prostheses are sterilized with gamma radiation in strictly controlled cycles. Open package directly before surgery to prevent contamination and observe the usual instructions for asepsis upon withdrawal.

Resterilization / Reprocessing

The prosthesis is intended for one-time use only. Reprocessing / resterilization is not permitted.

Storage

Storage conditions in unopened original packaging: Store in a dry place at room temperature, protected from direct sunlight. Brief fluctuations of temperature and humidity are permissible. Each prosthesis bears a batch number and an expiration date and may not be used after that date.

Documentation

It is recommended by the manufacturer to record the LOT-Number and Type of Prosthesis in the patient's medical record, operative record and implant pass by applying the adhesive labels.

Sales Restriction

US Federal Law restricts KURZ prostheses to sale by or on the order of a physician.

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EN

SYMBOLS

	Manufacturer		Item number
	Expiration date		Batch code
	Sterilized using irradiation		Do not use if package is damaged
	Keep away from sunlight		Do not resterilize
	Do not reuse!		Single sterile barrier system
	Store in a dry place		Single sterile barrier system with protective packaging inside
	Consult instructions for use		Single sterile barrier system with protective packaging outside
	Fragile, handle with care		Medical device
			Caution

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Heinz Kurz GmbH
KURZ Medical, Inc.