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STERILE R

CE 0124

IMPORTANT MEDICAL INFORMATION FOR BREATHE IMPLANT ÀWENGEN, IMPLANT FOR RHINOPLASTY

*Please read
carefully prior
to clinical
use!*

DESIGN	MATERIAL: TITANIUM	SIZE	mm	REF
	Breathe Implant àWengen	XS	3	6002 011
		S	4	6002 012
		M	5	6002 013
		L	6	6002 014
		XL	7	6002 015
		XXL	8	6002 016

INTENDED PURPOSE: For the expansion or stabilization of the nasal vestibule by means of spreading open the valvula (internal nasal valve) with the aim of achieving improved ventilation. **DESCRIPTION OF THE IMPLANT:** The product is a surgical invasive implant for long-term use. The implant is available in six sizes, it is made of pure titanium.

Indications: In cases of insufficient nasal breathing due to soft nasal stenosis resulting from anatomic conditions (nasal valve dysfunction) such as suction of the nasal sidewalls or inspiration or depression of the lower lateral nose.

Contraindications: Known allergies to titanium • Acute and chronic infectious diseases • General wound-healing impairments and bleeding tendency • Cases in which conservative treatment methods are sufficient

Intended User: The Breathe Implant à Wengen should only be used by surgeons with experience in Rhinology.

Intended Patient Group:

- All patients needing stabilization of the nasal vestibule by means of spreading the internal nasal valve are suitable for the operation after careful differential diagnostic measures.
- The implant can be used in all patients which have sizes of the nose in the range of the six implant sizes.
- As sufficient data for the treatment of children are not available, the implantation can only be recommended for patients older than 16 years.

Side Effects and Interactions: Interaction with microwaves or magnetic resonance imaging (MRI); see Warnings.

Possible complications:

- Postoperatively, there is a potential risk that the implant may extrude or that its contours may be overly visible, which would necessitate the removal of the implant. In such a case, the appropriate surgical measures must be initiated.
- Postoperatively, a thin layer of granulation tissue may form in the area immediately surrounding the implant during the healing phase.
- Postoperative infection is a complication of rhinoplasty and can occur early or late postoperatively.

Surgical Technique: Placement of the implant: The Breathe Implant must not extend below any part of the front edge of the lateral nasal cartilage, but must be placed approx. 1.5 mm above the edge. A cranial strip resection, as generally used for nasal tip reconstruction, must be avoided as this cartilage part of the lower lateral cartilage serves to partially cover and thus protect the Breathe Implant. The size and form of the lateral nasal cartilage are not modified, except if parts of the front edge, especially in the lateral area, press hard into the nasal vestibule. The valvula can be spread open according to the size of the selected implant. The Breathe Implant must not extend beyond the front edge of the cartilage. In a single-button seam technique, the Breathe Implant is sutured to the triangular cartilage through the front and rear slots of the implant using fine, non-absorbable sutures (5-0 or 6-0 Prolene or Nylon). The sutures must not be visible in the endonasal inspection.

Implant Selection: Intraoperative determination of the respective implant (XS - XXL) is possible with the aid of Sizers for Breathe Implants. The sizes XS, S, M, L, XL and XXL differ in the width of the median part. The intended purpose of the sizers is to have a tool to determine the correct size of the implant during the operation.

Accessories: Sizers for Breathe Implants: REF 8000249 (XS), 8000250 (S), 8000251 (M), 8000252 (L), 8000253 (XL), 8000254 (XXL)

Sterility: The Breathe Implants are provided sterile. Do NOT use, if the packaging is open or damaged prior to use! The Breathe Implants are sterilized with gamma radiation in strict accordance with the legal guidelines. The sterile packaging must only be opened immediately prior to use.

Packaging: Prosthesis packaging consists of:

- Plastic primary packaging (Polypropylene, PP) for the sterile implant.
- Sterile secondary packaging (GAG-PET and Tyvek®), inclusive 4 labels for patient's file, surgical record and implant card.
- Non-sterile storage packaging with instructions for use and implant card.

Resterilization/ Reprocessing:

- The Breathe Implants are single-use products. Resterilization/reprocessing is not permitted.
- For the sizers applies the "Reprocessing Instructions for Cleaning, Disinfection and Sterilisation of KURZ Instruments". The end of the product lifetime is significantly determined by the wear, the treatment process, the chemicals used and any possible damage caused by use. If disregarded, any liability is excluded.

MRI: The Breathe Implants are MR-Conditional for 1,5 T; 3,0 T and 7,0 T. For more detailed information about MRI please review www.kurzmed.com.

Warnings: The use of a damaged and / or soiled instrument (sizer for Breathe Implant) is the sole responsibility of the user. Placement of the implant: The Breathe Implant must not extend below any part of the front edge of the lateral nasal cartilage, but must be placed approx. 1.5 mm above the edge. The implant must not be bent.

The physician must inform the patient about following:

- Patients with metallic implants must not be exposed to microwave radiation;
- The patient should not pick their nose. Any annoying crusts should be removed by rinsing the nose with salt water and softening with nasal ointment;
- The patient should not pinch their nose tightly;
- In order to prevent the implant being bent, external forces e.g. participation in sports in which physical contact may occur intentionally or inadvertently must be avoided. Before wearing a respiratory protection mask or swim mask a physician must be consulted.
- In case of reddening, swelling or pain, an ENT specialist should be consulted to treat any potential soft tissue infections as soon as possible.

Storage: Store at dry room temperature in unopened original package. Each implant bears a batch number and an expiration date and must not be used after that date.

Disposal: Please refer to country-specific disposal guidelines.

Sales Restriction: The Breathe Implants are restricted to sale to or on the order of a specialist physician.

Documentation: The manufacturer recommends complete clinical, radiological and statistical documentation. Batch number, implant type and material must be entered into the patient's file, the surgical record or the implant card, using the adhesive labels provided with each implant.

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en

SYMBOLS

	Manufacturer		Item number
	Expiration date		Batch code
	Sterilized using irradiation		Do not use if package is damaged
	Keep away from sunlight		Do not resterilize
	Do not reuse!		Single sterile barrier system
	Store in a dry place		Single sterile barrier system with protective packaging inside
	Consult instructions for use		Single sterile barrier system with protective packaging outside
	Fragile, handle with care		Medical device
			Caution

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