NiTiBOND Stapes Prosthesis

(DE)	Wichtige Gebrauchsinformation	3
(EN)	Important Instructions for use.....	9
(FR)	Information importante avant utilisation	15
(IT)	Informazioni importanti per l'uso.....	21
(NL)	Belangrijke gebruiksaanwijzing.....	27
(PL)	Ważna informacja dla użytkownika.....	33
(PT)	Importante informação de utilização	39
(DK)	Vigtig brugsanvisning.....	45
(SE)	Viktig användarinformation	51
(ES)	Información importante para el uso.....	57
(TB)	Önemli tıbbi bilgiler.....	63
(CN)	重要的使用提示	69

WICHTIGE GEBRAUCHSINFORMATION

NiTiBOND STAPESPROTHESE

DE

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NiTiBOND Stapes Prosthesis

EN


Read carefully before clinical use!

Product Description

KURZ middle ear prostheses are intended for ossicular replacement to restore functionality to the middle ear in cases of pathological changes of the sound transmission system. The NiTiBOND stapes prosthesis acts as a bridge between the long process of the incus and the base of stapes, with the piston extending directly into the perilymphatic space of the inner ear. The loop consists of nickel-titanium alloy (Nitinol), while the piston and the shaft are made of pure titanium.


Attention!

The NiTiBOND stapes prosthesis is not suitable for placement with the piston on a section of tissue resulting in indirect extension into the perilymphatic space.


Piston-Diameter and Lengths

Ø 0.4mm		Ø 0.6mm	
REF	L (mm)	REF	L (mm)
1007 103	3.50	1007 153	3.50
1007 104	3.75	1007 154	3.75
1007 105	4.00	1007 155	4.00
1007 106	4.25	1007 156	4.25
1007 107	4.50	1007 157	4.50
1007 108	4.75	1007 158	4.75
1007 109	5.00	1007 159	5.00
1007 111	5.50	1007 161	5.50

Open Position

Closed Position

Intended Use

KURZ middle ear prostheses provide a partial or total replacement of the human ossicular chain when surgically indicated (see below). The objective is to transmit the greatest amount of sound energy possible to the inner ear without loss along the mechanical path.

Indications

Before a decision is made to implant a prosthesis, the patient's complete medical history must be reviewed:

- Otosclerosis (stapedial fixation) / congenital stapedial fixation
- Traumatic injury to the ossicular chain
- Malformation of the middle ear
- Revision surgery to correct inadequate hearing improvement, e.g. through dislocation of a prosthesis

Contraindications and Possible Risks

- Allergy to nickel, titanium, or nickel-titanium alloys (nitinol): A preoperative allergy test is recommended to minimize the risk of an allergic reaction. In patients with a known allergy to nickel, titanium, or nickel-titanium alloys (nitinol), do NOT use an implant containing these materials.
- Cases in which conservative treatment methods (e.g. a hearing device) are adequate.
- Acute middle ear inflammation which can lead to labyrinthitis or prosthesis dislocation.
- Inflammation of the auditory canal increases the risk of spreading infective foci into the middle or inner ear.
- If the ear undergoing surgery is the sole hearing source for the patient, there is a risk of complete bilateral deafness. In individual cases, this must be discussed between patient and physician.
- Acute and chronic infectious diseases
- General wound healing disorders

Possible Complications and Adverse Events

Adverse events or injuries can occur during or after the intervention. Extremely fine bony structures are manipulated and moved as part of the surgical procedure which can lead to iatrogenic trauma or an infection. These injuries may be irreversible or only correctable with revision surgery.

- Postoperative implant dislocation or fracture
- Tinnitus
- Necrosis or erosion in the vicinity of the contact zones (incus)
- Tissue irritation, scar formation, granulomas
- Recurring otitis media or progression of otosclerosis
- Dizziness
- Inner ear injury, including deafness
- Tympanic perforation
- Perilymphatic fistula
- Irritation or even damage to the facial nerve, including facial paralysis

Intolerance reactions such as granuloma formation and implant rejection, or interaction with medications are not known to date. In such cases, the prosthesis must be removed using conventional surgical procedures.

⚠ Warnings – Precautions

Intraoperative

Unintentional bending must be avoided during removal of the implant from the primary package; otherwise functional damage to the implant may occur. The middle ear structures must be manipulated as gently as possible. Prevention of injury to the inner ear must be ensured, especially in the region of the stapes and the oval window. Extreme care must be taken with the selection of the prosthesis length to reduce the risk of postoperative complications such as implant dislocation or the occurrence of inner ear symptoms, e.g. dizziness.

Postoperative

If a granuloma or perilymphatic fistula occurs postoperatively, initiate the required medical response immediately.


Measurement of Prosthesis Length

$$L_{\text{Prosthesis}} = L_1 + L_2$$

L1 = Measured distance

L2 = Immersion depth* (Piston)

*The determination of the immersion depth (L2) of the prosthesis is subject to the decision of the surgeon

Recommended Instrumentation

- KURZ measuring device (REF 8000 106) to determine the required length.

KURZ-Meter
(REF: 8000 106)

- Titanium forceps (REF 8000 136) for removing and handling the prosthesis.

Handling / Use / Implantation

Important: **Use only surgical laser for heat application.**

Figure 1. Overview – Contents


Figure 2 Remove Safety Pin

Push the Safety Pin downward with a pointed instrument until implant and Thermo-Dummy are released.


Figure 3 Raising the Thermo-Dummy

Raise the Thermo-Dummy by the eye with a small hook or forceps until it locks into upright position.


Figure 4 Laser Calibration

Calibrate the laser that will be used for heat application. The Thermo-Dummy is not an implant; it is intended exclusively for the determination of the exact power level required to close the NiTiBOND Stapes Prosthesis. Start with the lowest laser power setting and gradually increase to the setting that will close the Thermo-Dummy.


Figure 5 Prosthesis Removal

Remove the NiTiBOND Stapes Prosthesis with the titanium forceps. Make sure to grasp it at the midpoint of the piston to avoid prosthesis damage.



Figure 6 "Parking"

"Park" the implant in one of the customized slots within the primary box ($\varnothing$ 0.4 / 0.6 mm) to facilitate transfer to the middle ear.

Implantation

Figure 7 Transfer to the Middle Ear

Transfer from the primary package and insertion of the prosthesis into the middle ear should be performed, if at all possible, with a thin suction tube or fine micro forceps. Carefully apply the suction tube or micro forceps to the prosthesis loop.


Figure 8 Suspension on the Incus

Mount the prosthesis on the stapes foot plate by inserting the piston into the prepared opening/perforation. Then carefully slide the loop over the long incus process.



Figure 9 Laser Heat Application

Apply the laser heat sequentially to Points 1, 2 and 3 (red spots) to close the loop around the incus for an effective and safe implant connection.

It must be ensured that there is as little fluid as possible in the thermoactive zones as this could increase the amount of laser energy required.


Laser Procedure

Observe the sequence “Spot 1” – “Spot 2” – “Spot 3” during laser application! The setting should be equivalent to the result obtained during laser calibration (see Figure 4). Multiple laser shots may also be applied for each spot. However, due to different incus shapes, variations in the laser application sequence may be necessary. The laser heat setting for closing the loop starts with an As temperature ~55°C (131°F) and ends with Af ~ 70°C (158°F). The unique geometry of the thermal block assures that the heat does not reach the incus mucosa under normal circumstances.

Figure 10 Sealing the Stapes Prosthesis

If at all possible, insert the prosthesis upright into the perforation on the stapedial base plate. It is essential that the perforation is large enough to prevent loss of friction. Furthermore, the perforation must be sealed using tissue around the piston component. Once the prosthesis has been positioned and the perforation sealed, it is essential that the surgical coupling of the prosthesis to the incus is verified again. However, this may only be performed by way of indirect palpation of the incus or malleus, and not directly at the prosthesis loop.

Attention!

If implantation of the prosthesis is performed under local anesthesia (LA), the patient can be startled by the noise created during laser application. The patient should be forewarned of this to prevent unforeseeable head movements resulting in potential injury to the middle ear. If after closure of the loop rinsing is necessary for whatever reason, it should ONLY be done with solutions at body temperature.

If the procedure is performed with the aid of an endoscope, bear in mind that the use of an endoscope may result in temperatures that could cause the prosthesis loop to close prematurely.

MRI Safety Information

MR Conditional at 1.5 T; 3.0 T and 7.0 T, for detailed Information please refer to the MR Information on www.kurzmed.com

Warning:

The physician must inform the patient about the following: Patients with metallic implants are not to be subjected to microwave radiation. Severe variations in ambient pressure (scuba diving, diving head-first, explosions, fireworks, etc.) can damage the middle ear structures and/or cause hearing loss and/or equilibrium disorder and should be avoided.

Sterility

KURZ middle ear prostheses are provided sterile. Packaging contents are sterile, provided that the packaging was not damaged or opened. KURZ prostheses are sterilized with gamma adiation in strictly controlled cycles. Open package directly before surgery to prevent contamination and observe the usual instructions for asepsis upon withdrawal.

Resterilization/ Reprocessing

The prosthesis is intended for one-time use only. Reprocessing / resterilization is not permitted.

Storage

Storage conditions in unopened original packaging: Store in a dry place at room temperature, protected from direct sunlight. Brief fluctuations of temperature and humidity are permissible. Each prosthesis has a batch number and an expiration date. The prosthesis may not be implanted after the expiration date.

Documentation

It is recommended by the manufacturer to record the LOT-Number and Type of Prosthesis in the patient's medical record, operative record and implant pass by applying the adhesive labels.

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EN

SYMBOLS

	Manufacturer		Item number
	Expiration date		Batch code
	Sterilized using irradiation		Do not use if package is damaged
	Keep away from sunlight		Do not resterilize
	Do not reuse!		Single sterile barrier system
	Store in a dry place		Single sterile barrier system with protective packaging inside
	Consult instructions for use		Single sterile barrier system with protective packaging outside
	Fragile, handle with care		Medical device
			Caution

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**Heinz Kurz GmbH**
KURZ Medical, Inc.