AC Sizer® System Partial

REF 8000 540

(DE)	WICHTIGE GEBRAUCHSINFORMATION	2
(EN)	IMPORTANT INSTRUCTIONS FOR USE	4
(FR)	INFORMATIONS IMPORTANTES	6
(IT)	INFORMAZIONI IMPORTANTI PER L'USO	8
(ES)	INFORMACIÓN DE USO IMPORTANTE	10
(PT)	INFORMAÇÃO IMPORTANTE PARA USO	12
(RU)	ВАЖНЫЕ ИНСТРУКЦИИ ПО ИСПОЛЬЗОВАНИЮ	14
(CZ)	DŮLEŽITÉ INFORMACE O POUŽITÍ	16
(SRB)	VAŽNE KORISNIČKE INFORMACIJE	18
(PL)	WAŻNA INFORMACJA DLA UŻYTKOWNIKA	20
(IS)	MIKILVÆGAR UPPLÝSINGAR	22
(DK)	VIGTIG BRUGSANVISNING	24
(FI)	TÄRKEITÄ KÄYTTÖOHJEITA	26
(SE)	VIKTIG ANVÄNDARINFORMATION	28
(NO)	VIKTIG INFORMASJON OM BRUKEN	30
(TR)	KULLANIM İÇİN ÖNEMLİ TALİMATLAR	32
(NL)	BELANGRIJKE GEBRUIKSAANWIJZING	34
(GR)	ΣΗΜΑΝΤΙΚΕΣ ΠΛΗΡΟΦΟΡΙΕΣ ΓΙΑ ΤΟΝ ΧΡΗΣΤΗ	36
(HU)	FONTOS HASZNÁLATI INFORMÁCIÓK	38
(SI)	POMEMBNA NAVODILA ZA UPORABO	40
(LT)	SVARBI INFORMACIJA APIE NAUDOJIMĄ	42
(EE)	OLULINE KASUTUSINFORMATSIOON	44
(RO)	INFORMAȚII IMPORTANTE DE UTILIZARE	46
(CN)	重要的使用提示	48

STERILE R

CE₀₁₂₄

WICHTIGE GEBRAUCHSINFORMATION

DE AC Sizer® System Partial REF 8000 540

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IMPORTANT INSTRUCTIONS FOR USE

EN AC Sizer® System Partial REF 8000 540

Please read carefully before clinical application!

Intended purpose

The sizer disk is designed to allow the exact length of the middle ear prosthesis for implantation during a tympanoplasty, in which a partial ossicular chain reconstruction is performed, to be calculated. The plastic sizer disk is a sterile disposable product and must not be used more than once. Individual sizers are used for only a short time in the context of middle ear surgery and must not be left inside the middle ear for prolonged periods.

Using the sizer disk (partial)


fig. 1a


fig. 1b

1. Measurement

■ **Fig. 1a, 1b/2** The sizer disk for partial prostheses has six individual satellite-shaped sizers of different lengths (2.0-4.5 mm) (Fig. 1a, 1b). The relevant length is marked on the back of each sizer (Fig. 2). The cutouts and holes are designed to reduce the weight and to provide a clearer view when conducting measurements in the middle ear.

■ **Fig. 3** The sizer selected for the measurement process is cut out using microscissors and carried to the middle ear using an aspirator or other microsurgical instrument.

■ **Fig. 4** The bell-like part of the sizer is used to hold the head of the stapes during the measurement process. The plate-shaped part is designed to lie against the tympanic membrane. The thickness of the graft (e.g. disk of cartilage), later used to line the tympanic membrane is taken into account when calculating the length of the prosthesis. Once the optimum length has been determined, the partial prosthesis of the correct length can be selected / adjusted.

2. Expanding the bell on the prosthesis


■ **Fig. 5** If anatomical differences in the stapes require the bell on the partial prosthesis to be expanded, this can be achieved easily using the cone. To do so, press the bell carefully against the cone using a surgical needle or other microsurgical instrument. Under microscopic observation, the bell can be expanded accordingly.

3. Template for covering the prosthesis

■ **Fig. 1b/6** The integral templates on the back of the sizer disk (Fig. 1b) correspond to the minimum size of the graft for lining the tympanic membrane (Fig. 6).
Oval: for an oval prosthesis head plate
Round: for a round prosthesis head plate

4. Millimetre scale and sterile working plate

■ **Fig. 6** The millimetre scale is used for general comparative measurements, which might be useful during a middle ear operation. Basically, the back can also be used as storage space or a preparation field for various processes, such as working on the graft or for temporary storage.

Warnings / Precautions

During the surgical procedure, delicate bone structures are touched and moved, which can cause surgical trauma or infection. The middle ear structures must be treated as gently as possible. Trauma to the inner ear, e.g. by displacement of the stapes footplate, is particularly to be avoided in the region of the stapes and the oval window.

Sterility / Packaging

Packaging contents are sterile, provided that the packaging was not damaged or opened. This product is sterilised in strict compliance with the regulations governing gamma radiation. The sterile packaging should not be opened until immediately before use. When removing the product from the packaging, please observe the usual rules regarding asepsis. A storage pack contains 10 individual sterile packs.

Resterilization and reprocessing of AC Sizer® Systems

The AC Sizer® System is a disposable product. Reprocessing / resterilization is not permitted.

Storage:

Storage conditions in unopened original packaging: Store in a dry place at room temperature, protected from direct sunlight. Brief fluctuations of temperature and humidity are permissible. Do not use the product after the expiry date shown.

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.


fig. 2


fig. 3


fig. 4


fig. 5

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EN

SYMBOLS

	Manufacturer		Item number
	Expiration date		Batch code
	Sterilized using irradiation		Do not use if package is damaged
	Keep away from sunlight		Do not resterilize
	Do not reuse!		Single sterile barrier system
	Store in a dry place		Single sterile barrier system with protective packaging inside
	Consult instructions for use		Single sterile barrier system with protective packaging outside
	Fragile, handle with care		Medical device
			Caution

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**Heinz Kurz GmbH**
KURZ Medical, Inc.