## **AC Sizer® System Total**

### **REF 8000 550**

|       |                                             |    |
|-------|---------------------------------------------|----|
| (DE)  | WICHTIGE GEBRAUCHSINFORMATION .....         | 2  |
| (GB)  | (USA) IMPORTANT INSTRUCTIONS FOR USE .....  | 4  |
| (FR)  | INFORMATIONS IMPORTANTES .....              | 6  |
| (IT)  | INFORMAZIONI IMPORTANTI PER L'USO .....     | 8  |
| (ES)  | INFORMACIÓN DE USO IMPORTANTE .....         | 10 |
| (PT)  | INFORMAÇÃO IMPORTANTE PARA USO .....        | 12 |
| (RU)  | ВАЖНЫЕ ИНСТРУКЦИИ ПО ИСПОЛЬЗОВАНИЮ .....    | 14 |
| (CZ)  | DŮLEŽITÉ INFORMACE O POUŽITÍ .....          | 16 |
| (SRB) | VAŽNE KORISNIČKE INFORMACIJE .....          | 18 |
| (PL)  | WAŻNA INFORMACJA DLA UŻYTKOWNIKA .....      | 20 |
| (IS)  | MIKILVÆGAR UPPLÝSINGAR .....                | 22 |
| (DK)  | VIGTIG BRUGSANVISNING .....                 | 24 |
| (FI)  | TÄRKEITÄ KÄYTTÖOHJEITA .....                | 26 |
| (SE)  | VIKTIG ANVÄNDARINFORMATION .....            | 28 |
| (NO)  | VIKTIG INFORMASJON OM BRUKEN .....          | 30 |
| (TR)  | KULLANIM İÇİN ÖNEMLİ TALİMATLAR .....       | 32 |
| (NL)  | BELANGRIJKE GEBRUIKSAANWIJZING .....        | 34 |
| (GR)  | ΣΗΜΑΝΤΙΚΕΣ ΠΛΗΡΟΦΟΡΙΕΣ ΓΙΑ ΤΟΝ ΧΡΗΣΤΗ ..... | 36 |
| (HU)  | FONTOS HASZNÁLATI INFORMÁCIÓK .....         | 38 |
| (SI)  | POMEMBNA NAVODILA ZA UPORABO .....          | 40 |
| (LT)  | SVARBI INFORMACIJA APIE NAUDOJIMĄ .....     | 42 |
| (EE)  | OLULINE KASUTUSINFORMATSIOON .....          | 44 |
| (RO)  | INFORMAȚII IMPORTANTE DE UTILIZARE .....    | 46 |
| (CN)  | 重要的使用提示 .....                               | 48 |

**STERILE R**

**CE**<sub>0124</sub>

## WICHTIGE GEBRAUCHSINFORMATION

DE AC Sizer® System Total REF 8000 550

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## IMPORTANT INSTRUCTIONS FOR USE


**USA GB AC Sizer® System Total REF 8000 550**

**Please read carefully before clinical application!**

### Intended purpose

The sizer disk is designed to allow the exact length of the middle ear prosthesis for implantation during a tympanoplasty, in which a total ossicular chain reconstruction is performed, to be calculated. The plastic sizer disk is a sterile disposable product and must not be used more than once. Individual sizers are used for only a short time in the context of middle ear surgery and must not be left inside the middle ear for prolonged periods.

### Using the sizer disk (total)

Front

Back

#### 1. Measurement

■ **Fig. 1a, 1b/2** The sizer disk (Fig. 1a, 1b) for total prostheses has eight individual satellitishaped sizers of different lengths (3.0-7.0 mm). The relevant length is marked on the back of each sizer (Fig. 2). The cutouts and holes are designed to reduce the weight and to provide a clearer view when conducting measurements in the middle ear.

■ **Fig. 3** The sizer selected for the measurement process is cut out using microscissors and carried to the middle ear using an aspirator or other microsurgical instrument.

■ **Fig. 4** The stamp-like part of the sizer is used to support the stapes footplate during the measurement process. The plate-shaped part is designed to lie against the tympanic membrane. The thickness of the graft (e.g. disk of cartilage), later used to line the tympanic membrane is taken into account when calculating the length of the prosthesis. Once the optimum length has been determined, the total prosthesis of the correct length can be selected / adjusted.


Fig. 2


Fig. 3


Fig. 4

#### 2. Template for covering the prosthesis

■ **Fig. 1b/5** The integral templates on the back of the sizer disk (Fig. 1b) correspond to the minimum size of the graft for lining the tympanic membrane (Fig. 5). Oval: for an oval prosthesis head plate  
Round: for a round prosthesis head plate

#### 3. Millimetre scale and sterile working plate

■ **Fig. 5** The millimetre scale is used for general comparative measurements, which might be useful during a middle ear operation. Basically, the back can also be used as storage space or a preparation field for various processes, such as working on the graft or for temporary storage.


Fig. 5

### Warnings / Precautions

During the surgical procedure, delicate bones structures are touched and moved, which can cause surgical trauma or infection. The middle ear structures must be treated as gently as possible. Trauma to the inner ear, e.g. by displacement of the stapes footplate, is particularly to be avoided in the region of the stapes and the oval window.

### Sterility / Packaging

Packaging contents are sterile, provided that the packaging was not damaged or opened. This product is sterilised in strict compliance with the regulations governing gamma radiation. The sterile packaging should not be opened until immediately before use. When removing the product from the packaging, please observe the usual rules regarding asepsis. A storage pack contains 10 individual sterile packs.

### Resterilization and reprocessing of AC Sizer® Systems

The AC Sizer® System is a disposable product. Heinz Kurz GmbH Medizintechnik does not recommend the reprocessing and re-sterilization of used AC Sizer®. A product is considered 'used' when it has come into contact with foreign tissue or blood or when it has been removed from its packaging and touched, thereby undergoing mechanical stress. After cleaning, disinfection, resterilization and other handling, the AC Sizer® may potentially exhibit the following flaws:

- Insufficient sterility
- Insufficient non-pyrogenicity
- Presence of particles
- Formation of hazardous substances
- Microcracking
- Impairment of material properties

Due to the delicate designs and geometric structures of the AC Sizer®, the properties that guarantee product safety and effectiveness cannot easily be detected visually or without the appropriate testing equipment. Therefore, KURZ does not recommend reprocessing / resterilization of their devices. Should single-use prostheses undergo reprocessing for re-use, deterioration of the product quality must be anticipated. Liability for injury or damage caused by reprocessed implants/prostheses rests with the user and not with the manufacturer. KURZ cannot guarantee product safety and effectiveness after a single-use, CE-certified or FDA-cleared (or cleared / certified by another national regulatory authority) original product has been reprocessed / resterilized.

### Storage:

Storage conditions in unopened original packaging: Store in a dry place at room temperature, protected from direct sunlight. Brief fluctuations of temperature and humidity are permissible. Do not use the product after the expiry date shown.

**CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.**

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EN

## SYMBOLS

|                                                                                    |                                 |                                                                                     |                                                                       |
|------------------------------------------------------------------------------------|---------------------------------|-------------------------------------------------------------------------------------|-----------------------------------------------------------------------|
|  | Manufacturer                    |  | Item number                                                           |
|  | Expiration date                 |  | Batch code                                                            |
|  | Sterilized using irradiation    |  | Do not use if package<br>is damaged                                   |
|  | Keep away from sunlight         |  | Do not resterilize                                                    |
|  | Do not reuse!                   |  | Single sterile barrier<br>system                                      |
|  | Store in a dry place            |  | Single sterile barrier<br>system with protective<br>packaging inside  |
|  | Consult instructions<br>for use |  | Single sterile barrier<br>system with protective<br>packaging outside |
|  | Fragile, handle with care       |  | Medical device                                                        |
|                                                                                    |                                 |  | Caution                                                               |

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KURZ Medical, Inc.