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DESIGN	PRODUCT / TYPE	REFERENCE NUMBER	MATERIAL	WEIGHT (mg)	LENGTH (mm)	INDICATIONS FOR USE	IMPLANT AND INTENDED SITUATION
	ACP (Angular CliP Prosthesis)	1002 615 1002 617	Titanium (ASTM F67 Medical Grade)	4	2.25 3.25	For bridging defects at the long incudal process in patients with otherwise intact mobile ossicular chain.	The prosthesis consists of an angular stem with two, slightly staggered clamps for attachment to the incudal process and a clip that is 'clicked' onto the stapes head. No further manipulation is required.

Information for Use: Angular CliP Prosthesis (ACP)

Apply this Information for Use within the context of the "General Information for Ossiculoplasty Prostheses".

Surgical Technique:	The surgical techniques are demonstrated at various surgical seminars. Call the manufacturer/distributor for details or check programs of professional conferences. Attention: Removal from Package: Be absolutely sure to avoid accidental bending of implant during withdrawal from plastic container; otherwise proper function of the prostheses might be impeded.
Sterility:	KURZ middle ear prostheses are provided sterile. The prostheses are sterilized with gamma radiation in a strictly controlled cycle with sterility validation testing for each batch. Inspect package for punctures or other damage prior to surgery. Open storage package directly before surgery to prevent contamination.
Resterilization/ Reprocessing:	The prosthesis is intended for one-time use only. Resterilization / reprocessing is not permitted.
Warning:	The physician must inform the patient about the following: Patients with metallic implants are not to be subjected to microwave radiation. Severe variations in ambient pressure (scuba diving, diving head-first, explosions, fireworks, etc.) can damage the middle ear structures and/or cause hearing loss and/or equilibrium disorder and should be avoided.
MRI Safety Information:	MR Conditional at 1.5 T; 3.0 T and 7.0 T, for detailed Information please refer to the MR Information on www.kurzmed.com
Liability:	Complications or other adverse affects which might occur due to incorrect indication or unsatisfactory surgical technique, selection of unsuitable handling of them, asepsis, etc. come under the liability of the surgeon.
Storage:	Storage conditions in unopened original packaging: Store in a dry place at room temperature, protected from direct sunlight. Brief fluctuations of temperature and humidity are permissible. Each prosthesis bears a batch number and an expiration date and may not be used after that date.
Sale Restriction:	US federal law restricts KURZ prostheses to sale by or on the order of a physician.
Documentation:	It is recommended by the manufacturer to record the LOT-Number and Type of Prosthesis in the patient's medical record, operative record and implant pass by applying the adhesive labels.
Manufacturer:	Heinz Kurz GmbH, Tuebinger Strasse 3, 72144 Dusslingen, Germany Phone: +49 (0)7072 / 9179 - 0; Fax: +49 (0)7072 / 9179 - 79; Internet: www.kurzmed.com ; E-Mail: info@kurzmed.de USA: Kurz Medical Inc., 70 Chestnut Street Shrewsbury, MA 01545-USA Phone: 508.841.5900 Fax: 508.519.2672 Toll Free: 866.449.8020 Customer Service: orders@kurzmed.com

SYMBOLS

	Manufacturer		Item number
	Expiration date		Batch code
	Sterilized using irradiation		Do not use if package is damaged
	Keep away from sunlight		Do not resterilize
	Do not reuse!		Single sterile barrier system
	Store in a dry place		Single sterile barrier system with protective packaging inside
	Consult instructions for use		Single sterile barrier system with protective packaging outside
	Fragile, handle with care		Medical device
			Caution

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- Fabricante:** Heinz Kurz GmbH, Tübinger Straße 3, 72144 Dusslingen, Germany;
Wytwórca: Heinz Kurz GmbH, Tübinger Straße 3, 72144 Dusslingen, Germany,