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REF 1002 020

(DE)	WICHTIGE GEBRAUCHSINFORMATION.....	4
(EN)	IMPORTANT INSTRUCTIONS FOR USE	8
(FR)	INFORMATIONS IMPORTANTES.....	12
(IT)	INFORMAZIONI IMPORTANTI PER L'USO.....	16
(ES)	INFORMACIÓN DE USO IMPORTANTE	20
(PT)	INFORMAÇÃO IMPORTANTE PARA USO	24
(RU)	ВАЖНЫЕ ИНСТРУКЦИИ ПО ИСПОЛЬЗОВАНИЮ.....	28
(CZ)	DŮLEŽITÉ INFORMACE O POUŽITÍ	32
(SRB)	VAŽNE KORISNIČKE INFORMACIJE	36
(PL)	WAŻNA INFORMACJA DLA UŻYTKOWNIKA	40
(IS)	MIKILVÆGAR UPPLÝSINGAR.....	44
(DK)	VIGTIG BRUGSANVISNING	48
(FI)	TÄRKEITÄ KÄYTTÖOHJEITA	52
(SE)	VIKTIG ANVÄNDARINFORMATION	56
(NO)	VIKTIG INFORMASJON OM BRUKEN	60
(TR)	KULLANIM İÇİN ÖNEMLİ TALİMATLAR.....	64
(NL)	BELANGRIJKE GEBRUIKSAANWIJZING	68
(GR)	ΣΗΜΑΝΤΙΚΕΣ ΠΛΗΡΟΦΟΡΙΕΣ ΓΙΑ ΤΟΝ ΧΡΗΣΤΗ	72
(HU)	FONTOS HASZNÁLATI INFORMÁCIÓK	76
(SI)	POMEMBNA NAVODILA ZA UPORABO	80
(LT)	SVARBI INFORMACIJA APIE NAUDOJIMĄ	84
(EE)	OLULINE KASUTUSINFORMATSIOON	88
(RO)	INFORMAȚII IMPORTANTE DE UTILIZARE.....	92
(CN)	使用指南(临床使用前请仔细阅读)	96

STERILE R

CE₀₁₂₄

WICHTIGE GEBRAUCHSINFORMATION

DE TTP® VARIAC System Partial REF 1002 020

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IMPORTANT INSTRUCTIONS FOR USE

EN TTP® VARIAC System Partial REF 1002 020

Please read carefully before clinical application!

Description

The TTP® VARIAC System consists of two components: A TTP® VARIAC prosthesis made from titanium and a sizer disk made from plastic.

Sizer Disk
TTP® VARIAC Titanium Prosthesis in its
Plastic Box

Intended Use

The TTP® VARIAC titanium prosthesis is used in surgical indications as a partial replacement for the human ossicular chain in order to transmit sound energy mechanically from the tympanic membrane to the oval window of the cochlea with as little loss as possible. The sizer disk is designed to allow the exact length of the middle ear prosthesis for implantation during a tympanoplasty, in which a partial ossicular chain reconstruction is performed, to be calculated. The plastic sizer disk is a sterile disposable product and must not be used more than once. Individual sizers are used for only a short time in the context of middle ear surgery and must not be left inside the middle ear for prolonged periods.

Indication

Before deciding to implant a passive middle ear prosthesis to reconstruct the human ossicular chain, it must be ensured that on the one hand the complete history of the patient is also taken into account, and on the other hand that the ear microscopy findings as well as the intra-operative situation after opening up the middle ear and exposing the damaged or altered structures suggest that use of a prosthesis is indicated. This decision can, according to the present-day state of medical knowledge and technology, only be reached during ongoing surgery after the operative site has been exposed.

Implantation of a passive middle ear prosthesis to reconstruct the human ossicular chain can be indicated on account of the following underlying diseases:

- Chronic middle ear inflammation (also following removal of a tumour, e.g. cholesteatoma) with functional impairment of the ossicular chain
- Traumatic injury to the ossicular chain
- Malformation of the middle ear
- Revision surgery due to inadequate hearing improvement, e.g. due to dislocation of a (previously implanted) prosthesis

Contraindications

- Cases in which conservative treatment methods are adequate
- Acute otitis media, which can lead to dislocation of the prosthesis
- Acute and chronic infectious diseases
- Titanium allergy
- General wound-healing disorders

Possible complications / Unanticipated adverse effects

Possible lasting consequences / damage may occur during or after the operation. During the surgical procedure, delicate bone structures are touched and moved, which can cause surgical trauma or infection. This damage may be irreversible, or correctable only by a revision operation.

- Postoperative dislocation of the prosthesis
- Recurrent otitis media
- Giddiness
- Tissue irritation / scar formation / granuloma
- Perforation of the eardrum

- Adhesions around the prosthesis
- Various degrees of inner ear injury up to deafness
- Tinnitus
- Irritation or even damage to the facial nerve, including facial paralysis

Warnings / Precautions

Intraoperative

When removing the product from the primary packaging, please avoid inadvertently bending it, as this may cause functional damage to the implant. The middle ear structures must be treated as gently as possible. Trauma to the inner ear, e.g. by displacement of the stapes footplate, is particularly to be avoided in the region of the stapes and the oval window. The length of the prosthesis to be implanted should be selected with the utmost attentiveness in order to avoid potential problems, e.g. subsequent dislocation of the implant. The prosthesis head plate must be covered by a 0.5 mm thick graft (e.g. disk of cartilage) between it and the tympanic membrane. When selecting the length of prosthesis, please read the sizer disk instructions (see below).

Postoperative

If granulomas form or an allergic reaction occurs postoperatively, the necessary medical steps must be implemented immediately. The patient should be instructed in hygienic after-care of the operated ear.

Safety Information

1. MR Conditional at 1.5 T; 3.0 T and 7.0 T, for detailed Information please refer to the MR Information on www.kurzmed.com

Note

The doctor should inform the patient of the following: Patients with metal implants must not be irradiated with microwaves. Severe variations in ambient pressure (scuba diving, diving head-first, explosions, fireworks, etc.) can damage the middle ear structures and / or cause hearing loss and / or equilibrium disorder and should be avoided.

Using the sizer disk


Fig. 1a
Front


Fig. 1b
Back

1. Measurement

■ **Fig. 1a, 1b/2** The sizer disk for partial prostheses has six individual satelliteshaped sizers of different lengths (2.0-4.5 mm) (Fig. 1a, 1b). The relevant length is marked on the back of each sizer (Fig. 2). The cutouts and holes are designed to reduce the weight and to provide a clearer view when conducting measurements in the middle ear.


Fig. 2

■ **Fig. 3** The sizer selected for the measurement process is cut out using microscissors and carried to the middle ear using an aspirator or other microsurgical instrument.




Fig. 3

■ **Fig. 4** The bell-like part of the sizer is used to hold the head of the stapes during the measurement process. The plate-shaped part is designed to lie against the tympanic membrane. The thickness of the graft (e.g. disk of cartilage), later used to line the tympanic membrane is taken into account when calculating the length of the prosthesis. Once the optimum length has been determined, the partial prosthesis of the correct length can be selected.


2. Adjusting the length on the prosthesis


■ **Fig. 5** The TTP® VARIAC prosthesis is inserted in the purpose-designed recess in the sizer disk according to the length calculated beforehand (see below: Using the sizer disk). The slots on the bell of the prosthesis are aligned in such a way that the bell slides downwards along the bridging pins. This procedure should be carried out using the titanium tweezers (REF 8000 136) supplied separately.



■ **Fig. 6** The bell of the prosthesis must be positioned at the base of the recess and the head plate must be seated fully in the purpose-designed cutout so that the upper surface of the head plate engages with the upper surface of the sizer disk.

■ **Fig. 7/8** To close that head plate the Titanium Micro Closing Forceps (REF 8000 137) must be used. The handle inside the head plate of the prosthesis is straightened using the closing forceps (REF 8000 137), fixing the head plate to the shaft of the prosthesis (Fig. 7). When doing this, care should be taken to ensure that the correct side of the closing forceps is in the corresponding position on the head plate of the prosthesis. "INSIDE" for the inside of the head plate and "OUTSIDE" for the outside of the head plate (Fig. 8a, 8b).


■ **Fig. 9** Finally, the shaft is cut off next to the head plate using the cutting forceps (REF 8000 171). The pin that this leaves on the head plate is later used to support the fixation of the graft between the head plate and the tympanic membrane, and prevents postoperative dislocation. The length of the pin can be varied, but it should at any rate correlate with the thickness of the intended graft.


3. Expanding the bell on the prosthesis

■ **Fig. 10** If anatomical differences in the stapes require the bell on the partial prosthesis to be expanded, this can be achieved easily using the cone. To do so, press the bell carefully against the cone using a surgical needle or other microsurgical instrument. Under microscopic observation, the bell can be expanded accordingly.


4. Template for covering the prosthesis

■ **Fig. 1b/11** The integral templates on the back of the sizer disk (Fig. 1b) correspond to the minimum size of the graft for lining the tympanic membrane (Fig. 11).
Oval: for an oval prosthesis head plate
Round: for a round prosthesis head plate

5. Millimetre scale and sterile working plate

■ **Fig. 11** The millimetre scale is used for general comparative measurements, which might be useful during a middle ear operation. Basically, the back can also be used as storage space or a preparation field for various processes, such as working on the graft or for temporary storage.

Sterility

KURZ middle ear prostheses are provided sterile.

Packaging contents are sterile, provided that the packaging was not damaged or opened. This product is sterilised in strict compliance with the regulations governing gamma radiation. The sterile packaging should not be opened until immediately before use.

Resterilization and reprocessing

The TTP-VARIAC System is intended for one-time use only. Reprocessing / resterilization is not permitted.

Storage:

Storage conditions in unopened original packaging: Store in a dry place at room temperature, protected from direct sunlight. Brief fluctuations of temperature and humidity are permissible. Each prosthesis bears a batch number and an expiration date and may not be used after that date.

Documentation

It is recommended by the manufacturer to record the LOT-Number and Type of Prosthesis in the patient's medical record, operative record and implant pass by applying the adhesive labels.

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

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EN

SYMBOLS

	Manufacturer		Item number
	Expiration date		Batch code
	Sterilized using irradiation		Do not use if package is damaged
	Keep away from sunlight		Do not resterilize
	Do not reuse!		Single sterile barrier system
	Store in a dry place		Single sterile barrier system with protective packaging inside
	Consult instructions for use		Single sterile barrier system with protective packaging outside
	Fragile, handle with care		Medical device
			Caution

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Heinz Kurz GmbH
KURZ Medical, Inc.