


CARTILAGE PUNCH

- » Punches out a cartilage shoe for stabilization of a KURZ total prosthesis on the stapes footplate
- » Centers position of a prosthesis to increase acoustical transmission
- » Easy and fast to reduce surgical time



MIDDLE EAR INTELLIGENCE





Total Prosthesis in-situ

Image courtesy of K. B. Hüttenbrink,
Cologne, Germany

*The Cartilage Punch cannot be used together with the
Regensburg Type Total Prosthesis*

ITEM	REF
Cartilage Punch incl. Instrument Tray	8000 200
REPLACEMENT PARTS	REF
Instrument Tray incl. Inlay	8000 176
Spring	8000 198
Plastic Base POM (black)	8000 207

Material: Sterilizable instrument quality stainless steel

Heinz Kurz GmbH Medizintechnik | Tuebinger Strasse 3
72144 Dusslingen | Germany | Phone: +49 (0)7072/9179-0
E-Mail: info@kurzmed.de | www.kurzmed.com

KURZ Medical, Inc. | 5126 South Royal Atlanta Drive Tucker,
GA 30084 | USA Toll Free Order #: 866.449.8020 | Phone:
770.349.6330 | E-Mail: info@kurzmed.com | www.kurzmed.com



02/2018-M9600374

Prescription Devices

Caution: United States Federal Law restricts devices identified herein to sale by or on the order of a physician. This brochure does not replace nor does it set forth the complete contents of the "Instructions for Use" for the products in this brochure, and is not a substitute for reviewing and understanding that important information. Therefore, before using any products included in this brochure, please review the entire contents of the respective "Instructions for Use" information.

Because of registration the products are currently not available in all countries.



MIDDLE EAR INTELLIGENCE