

KURZ PRECISE CARTILAGE KNIFE

- » Precise cartilage slices from 0.1-0.7 mm
- » Multiple usage in middle ear surgery
- » Easy handling reduces surgical time



MIDDLE EAR INTELLIGENCE





ITEM	REF
KURZ Precise Cartilage Knife Set (blade holder, cutting block, distance plates, blades) incl. Instrument Tray	8000 155
Distance Plate for Precise 1mm (for compression of fascia)	8000 105

Supplies	REF
Blades sterile single packed, 10 pieces	8000 140

Material:
Sterilizable instrument quality stainless steel

REPLACEMENT PARTS	REF
1 Set of Distance Plates (0.1/0.2/0.3 mm)	8000 102
KURZ Precise Blade Holder incl. Knurled Screw M3, St. 1.4305	8000 103
Instrument Tray KURZ Precise	8000 177
Knurled Screw M3, St. 1.4305	8000 190
Nut Precise Cutting Block M6	8000 191
KURZ Precise Cutting Block incl. Nut Precise Cutting Block M6	8000 110
Instrument Tray KURZ Precise incl. Plate POM+Tray stainless steel	8000 144



02/2018-M9600373

Prescription Devices

Caution: United States Federal Law restricts devices identified herein to sale by or on the order of a physician. This brochure does not replace nor does it set forth the complete contents of the "Instructions for Use" for the products in this brochure, and is not a substitute for reviewing and understanding that important information. Therefore, before using any products included in this brochure, please review the entire contents of the respective "Instructions for Use" information.

Because of registration the products are currently not available in all countries.

Heinz Kurz GmbH Medizintechnik | Tuebinger Strasse 3
72144 Dusslingen | Germany | Phone: +49 (0)7072/9179-0
E-Mail: info@kurzmed.de | www.kurzmed.com

KURZ Medical, Inc. | 5126 South Royal Atlanta Drive Tucker,
GA 30084 | USA Toll Free Order #: 866.449.8020 | Phone:
770.349.6330 | E-Mail: info@kurzmed.com | www.kurzmed.com



MIDDLE EAR INTELLIGENCE